

# VASYLI<sup>®</sup> MEDICAL

ENHANCED BIOMECHANICAL FUNCTION

  
PATTERSON  
MEDICAL



HELP RESTORE NATURAL FOOT FUNCTION WITH  
VASYLI ORTHOTICS & VIONIC/ORTHAHEEL FOOTWEAR

Patterson Medical  
**EXCLUSIVE**


*Performance through innovation*

**VASYLI<sup>®</sup>**

**Footwear & orthotics help restore natural foot function and well-being with every step.**

Our feet were designed to walk on soft, natural surfaces like soil and sand. Yet today we spend most of our time standing, walking and running on unnatural, hard surfaces.

#### Natural relief



Before (left) After (right)

Clinical research supports Vasyli's success in restoring the foot's natural (neutral) alignment, enabling the feet, ankles and legs to function dynamically as they were intended. Vasyli's unique orthotic is built in to our collection of footwear.

## Orthotics - A Natural Cure...

It is estimated that the majority of the population over pronates due to the hard, flat unnatural surfaces we walk on daily.

Over pronation often disrupts normal knee function and hip alignment and increases forces on the lower back muscles.

These biomechanical deficiencies are common in all types of patients from children to the elderly, from top athletes to people with a sedentary lifestyle.

By realigning the lower limb to its natural angle Vasyli® orthotics ensure correct foot function and help relieve common biomechanical complaints.

In combination with the practitioner's regular treatment regime, orthotics are usually successful in alleviating chronic pains.



the Problem...



the Cure

#### Ideal for:

- Plantar Fasciitis
- Metatarsalgia
- Tired/Aching Legs
- Lower Back Pain
- Ilio-tibial Band Syndrome
- Achilles Tendonitis
- Hallux Abducto Valgus
- Patello Femoral Pain
- Tibial Stress Syndrome

“Supported by science, proven by customer satisfaction.”

The Vasyli Think Tank™ is at the very heart of Vasyli.

## Vasyli Custom Easyfit Ideal for Multiple Shoe Styles

### Blue = Medium Density


This model features a lateral cut-away providing an easy fit into hard-to-fit footwear, including men's moccasins, loafers, and ladies (low-heeled) court shoes. Sold as pairs. Latex free. See chart below for ordering information.



**Phillip J. Vasyli**  
World renowned podiatrist, Chairman and founder of Vasyli International has always been committed to delivering quality treatments and education designed to enhance practitioners' knowledge and patient's health and well-being.

## Vasyli Dananberg Orthotic Treating 1st Ray Function Just Got Easier!

Developed with podiatrist Howard Dananberg, renowned as a world leader in concepts relating to 1st ray function, the Vasyli Dananberg orthotic is a fast and simple alternative to a 1st Ray cutout. The 2 removable plugs allow you to customize the orthotic to assist in maximizing 1st ray function. In addition, this orthotic is also heat moldable to personalize the rearfoot control. Sold as pairs. Latex free. See chart below for ordering information.



Sorbon shock dot designed to dissipate force along the angle of strike.

2 removable plugs allowing you to customize the orthotic further to assist in maximizing first ray function.

Heat-moldable EVA allowing customization of the orthotic.




**Howard J. Dananberg**  
Podiatrist & member of the Vasyli Think Tank

## The Vasyli McConnell Slimfit Orthotic has just made fitting orthotics into women's high heels and sandals much easier

Developed with physical therapist Jenny McConnell, renowned as a world leader in treatment of the lower limb and taping techniques. With a new slimmer profile and enhanced metatarsal support, the new Vasyli McConnell slimfit device "disappears" under the foot making it the perfect solution for hard-to-fit footwear.

- Enhanced metatarsal control
- Slim design fits most heeled shoes & sandals
- Dual-density compound
- Supports & realigns the foot

Easyfit	Dananberg	McConnell	Size	Men's Shoe Size	Women's Shoe Size
—	—	—	Kids Sm	Youth 12-13	Youth 12-13
5609-31-01	—	—	Kids Lg	1½ - 3	1½ - 4
5609-31-02	5630-11-01	5629-89-01	XS	3½ - 5	4½ - 6
5609-31-03	5630-11-02	5629-89-02	S	5½ - 7	6½ - 8
5609-31-04	5630-11-03	5629-89-03	M	7½ - 9	8½ - 10
5609-31-05	5630-11-04	5629-89-04	L	9½ - 11	10½ - 12
5609-31-06	5630-11-05	—	XL	11½ - 13	—



**Jenny McConnell**  
Physiotherapist & member of the Vasyli Think Tank

The Vasyli Think Tank™ is at the very heart of Vasyli.



## Brian Hoke

Received his degree in Physical Therapy from Indiana University and earned his post-professional Doctor of Physical Therapy from Boston University. He is board certified as a Sports Clinical Specialist by the American Board of Physical Therapy Specialties. Brian is co-owner and Director of Atlantic Physical Therapy, a physical therapist owned private practice in Virginia Beach, Virginia. Brian lectures extensively in the USA and internationally.



## Vasyli® Hoke Supination Control Orthotic

By utilizing a triplanar (three dimensional) perspective to normalize foot biomechanics, this orthotic effectively addresses forefoot valgus, plantarflexed 1st ray and forefoot equinus.

The excessive supinator typically has:

- In the Sagittal Plane, a forefoot equinus frequently accompanied by a plantarflexed 1st ray
- In the Frontal Plane, a forefoot valgus
- In the Transverse Plane, a forefoot adductus with a prominent fifth metatarsal base.


The Vasyli Hoke device integrates triplanar features to address each of these anatomical considerations.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0816-09-957	XS	3½ – 5	4½ – 6
0816-09-965	S	5½ – 7	6½ – 8
0816-09-973	M	7½ – 9	8½ – 10
0816-09-981	L	9½ – 11	10½ – 12
0816-09-999	XL	11½ – 13	N/A



## Thomas G. McPoil

Currently a Professor in the School of Physical Therapy at Regis University in Denver, CO. He has specialized in the evaluation and management of foot and ankle disorders since 1979. His scholarly contributions have systematically examined foot and ankle function from both scientific and clinical perspectives.



## Vasyli® McPoil Tissue Stress Relief Orthotic

During activity, a certain degree of foot pronation is necessary for normal lower extremity mobility. To accommodate for this individual specific degree of "normal" foot pronation and simultaneously prevent over pronation, the Vasyli McPoil Tissue Stress Relief Orthotic incorporates a dual density design that achieves both functional control and excellent compliance. The dual-density, total arch contact principle of the Vasyli McPoil foot orthotic design is applicable in the treatment regimes of both athletes and non-athletes with lower extremity conditions such as: anterior knee pain, medial tibial stress syndrome, posterior tibialis tendinopathy, chronic plantar heel pain ("plantar fasciitis/fasciosis"), general arch pain, and mid-foot pain.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0816-10-005	XS	3½ – 5	4½ – 6
0816-10-013	S	5½ – 7	6½ – 8
0816-10-021	M	7½ – 9	8½ – 10
0816-10-039	L	9½ – 11	10½ – 12
0816-10-047	XL	11½ – 13	N/A



Influencing product development to benefit both practitioner and patient alike.

**Wider and thinner fore foot profile than the original Armstrong Orthotic!**



Top view

Bottom view


## Vasyli Armstrong II Sensitive Feet Orthotic

Whether you are dealing with an insensate foot, a hyper-sensitive foot type, or a diabetic foot, reducing pressure and shear forces should be playing a significant roll in your patients orthotic treatment regime.

Developed with podiatrist David Armstrong, renowned as a world leader in diabetic foot care, the Vasyli Armstrong II is the most advanced foot orthotic in its class. Featuring patented Glidesoft technology; the dorsal surface of the Vasyli-Armstrong II orthotic moves with the plantar surface of the foot dramatically reducing shear/frictional forces. This in-turn assists in the reduction of skin lesions and the potential for infection/ulceration.

The VAI is re-inventing the way leading health professionals care for sensitive feet! Ideal for the diabetic foot, patients suffering from rheumatoid arthritis or plantar fat pad atrophy causing painful sensitivity issues. This orthotic offers a clinically validated and scientifically tested shear and pressure reduction solution. Customize by using a standard heat gun and/or adding Vasyli postings. Sold as pairs. Latex free.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0815-80-406	XS	4-5	5-6
0815-80-414	S	6-7	7-8
0815-80-422	M	8-9	9-10
0815-80-430	L	10-11	11-12
0815-80-448	XL	12-13	13-14
0815-80-455	XXL	14-15	—



**David Armstrong**  
Professor of Surgery  
Director of Lower  
Extremity Research and  
the Southern Arizona Limb  
Salvage Alliance (SALSA)  
University of Arizona,  
Tucson. Member of the  
Vasyli Think Tank.



Top view

Bottom view

- The heat moldable high density shell (A) extends through the lateral aspect of the fore foot assisting in locking up the mid tarsal joint enhancing the propulsive phase of gait
- The PU fore foot (B) provides both comfort and function
- The 1st ray plantarflex zone (C) facilitates fore foot stability and Hallux dorsiflexion
- The PU shock dot in the heel (D) assists in attenuating heel strike forces

## Vasyli Prior Sports Orthotic

The Vasyli + Prior incorporates a subtle 4mm heel raise—perfect for helping with tight posterior muscles and associated conditions that are often exacerbated in flat lasted sports footwear. The rear foot post of 3° + the 4mm raise forms a smart medial wedge for control without the “bulk.” Low profile, medial arch contour ensures comfort; the high density, heat moldable shell provides further customization, enhancing fit and performance.

The Vasyli Prior Sports Orthotic embodies the essential elements of function and control while satisfying the fit characteristics required in low profile / hard to fit performance footwear. Basically any boot of shoe with a stud, cleat or spike.

Item No.	Size	Men's Shoe Size	Women's Shoe Size
0815-80-463	XS	3½ – 5	4½ – 6
0815-80-471	S	5½ – 7	6½ – 8
0815-80-489	M	7½ – 9	8½ – 10
0815-80-497	L	9½ – 11	10½ – 12
0815-80-505	XL	11½ – 13	—



**Trevor Prior**  
Director of Premier podiatry  
Limited and P2L, Trevor currently  
consults as a Podiatric  
Surgeon at Homerton University  
Hospital in the UK. Member  
of the Vasyli Think Tank.



## Vasyli Shock Absorber

Has a mild intrinsically posted rearfoot to help reduce excess pronation and is made from soft PU and shock absorbing materials to cushion the foot. Whilst not as corrective as the standard Vasyli custom orthotics it does offer mild control not normally found in shock absorbing insoles. Built-in metatarsal raise, supporting the forefoot area. Great for supinators and can be posted using any of the Vasyli postings and additions. As a result, it gives your patient excellent comfort, biomechanical control and lasting pain relief. Sold as pairs. Latex free.

**Full Length Only**



## Vasyli Custom Extended Slimfit

Specifically designed to fit ladies shoes, such as fashion shoes and high-heels. The flexible base allows easy contouring of the shank of the shoe. The rearfoot angle is reduced to 2 degrees due to the high heel factor (tibio talar axis) whose action inverts the foot as the heel is raised. A built-in metatarsal raise ensures correct alignment of the metatarsal bones, helping to prevent Metatarsalgia and Morton's Neuroma. Sold as pairs. Latex free. See chart for ordering information.

**Sulcus Length Only**



## Green = Low Density

### Available in 3/4 Length & Full Length

Ideal for elderly, pediatric and arthritic patients who require special consideration. Patients adapt easily to the mild functional control while the correct biomechanical position is maintained. Sold as pairs. Latex free.



## Blue = Medium Density

### Available in 3/4 Length & Full Length

Designed for the "average" patient with general orthotic needs, blue density orthotics provide excellent biomechanical control. Sold as pairs. Latex free.



## Red = High Density

### Available in 3/4 Length & Full Length

Designed for high performance, this range provides both durability and solid biomechanical control. Red, high performance orthotics are ideal for athletes, high activity users and overweight patients. Sold as pairs. Latex free.

Green		Blue		Red		Shock Absorbers	Extended Slimfit	Size	Men's Shoe Size	Women's Shoe Size
3/4th Length	Full Length	3/4th Length	Full Length	3/4th Length	Full Length					
5609-34-01	5609-33-01	5609-29-01	5609-27-01	5609-30-01	5609-28-01	—	—	KS	Youth 12–13	Youth 12–13
5609-34-02	5609-33-02	5609-29-02	5609-27-02	5609-30-02	5609-28-02	—	—	KL	1½ – 3	1½ – 4
5609-34-03	5609-33-03	5609-29-03	5609-27-03	5609-30-03	5609-28-03	5629-94-05	5609-32-01	XS	3½ – 5	4½ – 6
5609-34-04	5609-33-04	5609-29-04	5609-27-04	5609-30-04	5609-28-04	5629-94-01	5609-32-02	S	5½ – 7	6½ – 8
5609-34-05	5609-33-05	5609-29-05	5609-27-05	5609-30-05	5609-28-05	5629-94-02	5609-32-03	M	7½ – 9	8½ – 10
5609-34-06	5609-33-06	5609-29-06	5609-27-06	5609-30-06	5609-28-06	5629-94-03	5609-32-04	L	9½ – 11	10½ – 12
5609-34-07	5609-33-07	5609-29-07	5609-27-07	5609-30-07	5609-28-07	5629-94-04	—	XL	11½ – 13	—
5609-34-08	5609-33-08	5609-29-08	5609-27-08	5609-30-08	5609-28-08	—	—	XXL	13½ – 15	—

## Vasyli Heel Pain Orthotic

Specifically designed to provide instant relief from painful heel conditions. The Heel Piece is also designed for hard-to-fit footwear such as ladies fashion shoes and men's and women's dress shoes where other orthotics may be too bulky. Sold as pairs. Latex free.

**5609-67-01** Small, Fits shoe size: Men's 5-8; Women's 6-9

**5609-67-02** Large, Fits shoe size: Men's 9-13; Women's 10+



## Vasyli Orthotic Wedges and Additions

A full range of self-adhesive VASYLI orthotic wedges and additions have been designed to assist in creating a better position for the foot structure. The fact is all feet are not the same and often there's a need to modify the orthotic for specific problems that require extra biomechanical control. Package of 10. Latex free.


**5609-39** Metatarsal Raise, pkg. of 10

Wedges			Degree
Forefoot		Rearfoot	
Valgus	Varus	Varus	
—	<b>5609-37-01</b>	<b>5609-35-01</b>	2
<b>5609-36-01</b>	<b>5609-37-02</b>	<b>5609-35-02</b>	4
<b>5609-36-02</b>	<b>5609-37-03</b>	—	6

Heel Lift			Lift
Small	Medium	Large	
<b>5609-38-01</b>	<b>5609-38-04</b>	<b>5609-38-07</b>	4 mm
<b>5609-38-02</b>	<b>5609-38-05</b>	<b>5609-38-08</b>	6 mm
<b>5609-38-03</b>	<b>5609-38-06</b>	<b>5609-38-09</b>	8 mm

**Mixed Addition Pack.** This set provides you with a sampling of the Vasyli additions. Includes 1 pair of heel lifts, forefoot varus, rearfoot varus, forefoot valgus and metatarsal raisers. † **5629-95**

†Size and degrees vary per pack.



## Vasyli® Medical Starter Packs.

Purchase everything you need to get started making custom insoles, for your patient population, in one convenient kit. **All kits include a free heat gun and manual.**



### Active Starter Pack.

Includes 4 pairs of Dananberg insoles (sizes xsmall through large); 3 pairs full-length blue insoles (sizes small, medium and large); 2 pairs full-length red insoles (sizes large and xlarge); 1 pkg. Rearfoot Varus 2 degree (item #56093501) and 2 packs of mixed additions pkg (item #562995).

### All Rounder Starter Pack.

Includes 6 pairs of Dananberg insoles (2 pr of each sizes small through large); 3 pairs 3/4-length blue insoles (sizes small, medium and large); 2 pairs Extended SlimFit insoles (sizes xsmall and small); and 2 packs of mixed additions pkg (item #562995).

### The Sports Starter Pack.

Includes 5 pairs of Dananberg insoles (sizes xsmall through xlarge); 5 pairs of Prior insoles (sizes xsmall through xlarge) and 2 packs of mixed additions pkg (item #562995).

**0816-12-167** Active Starter Pack

**0816-12-175** All Rounder Starter Pack

**0816-12-183** Sports Starter Pack

### Tools

**Vasyli Heat Gun.**  
**0815-23-547**



**0815-23-570**

**Economy Magnetic Inclinator.**

Same tool that is used in Vasyli workshops is now available for your convenience.

**0815-23-570**



# Vionic Sandals for Men

The ideal motion control footwear for after-sport recovery.

**Vionic**  
with orthaheel™



Durable VIBRAM outsole for excellent traction on various surfaces. Biomechanically contoured, motion-control footbed made from triple-density EVA for stability in the rear foot and cushioning in the heel and forefoot. EVA midsole and VIBRAM outsole. Received APMA Seal of Acceptance.



## Vionic Bryce. ⚠

Performance men's toe post sandal with flexible leather uppers, padded jersey liner and a soft toe post for a comfortable fit. Available in whole sizes only, medium (D) width. Latex free.

Black/Grey	Khaki	Men's Shoe Size
0815-89-324	0815-89-399	7
0815-89-332	0815-89-407	8
0815-89-340	0815-89-415	9
0815-89-357	0815-89-423	10
0815-89-365	0815-89-431	11
0815-89-373	0815-89-449	12
0815-89-381	0815-89-456	13



## Vionic Boyes. ⚠

Performance men's adjustable back-strap sandal offers secure comfort for active lifestyles. Flexible leather uppers feature fully adjustable straps with hook-and-loop closure and soft padded jersey liner. Available in whole sizes only, medium (D) width. Latex free

Black	Taupe	Men's Shoe Size
0815-96-535	0815-96-592	8
0815-96-543	0815-96-600	9
0815-96-550	0815-96-618	10
0815-96-568	0815-96-626	11
0815-96-576	0815-96-634	12
0815-96-584	0815-96-642	13



## Vionic Base. ⚠

Provides the foundation for stable sports recovery. The slip-on, toe post sandal is built upon a triple-dense midsole, and features a pop of color etched in architectural detailing. Available in whole sizes only, medium (D) width. Latex free.

Black/Silver	Choc/Red	Men's Shoe Size
0815-96-329	0815-96-386	8
0815-96-337	0815-96-394	9
0815-96-345	0815-96-402	10
0815-96-352	0815-96-410	11
0815-96-360	0815-96-428	12
0815-96-378	0815-96-436	13

⚠ This product contains natural rubber latex or dry natural rubber.



## Vionic Traverse. ⚠

This slide-on, adjustable sandal crosses conventional boundaries to provide the comfort, support and stability demanded by athletes of all levels. Available in whole sizes only, medium (D) width. Black/Silver only. Latex free.

Item No.	Men's Shoe Size
0815-96-444	8
0815-96-451	9
0815-96-469	10
0815-96-477	11
0815-96-485	12
0815-96-493	13

Half sizes for men should go up a size (i.e. 10½ goes to 11)  
Half sizes for women should go down a size (i.e. 7½ goes to 7)



# Vionic Sandals for Women

The ideal motion control footwear for after-sport recovery.



## Vionic Cascade. ⚠

Performance women's toe post sandal with soft leather uppers, padded jersey liner with bungee strap overlay and a soft nylon webbing toe post. Available in whole sizes only in medium/wide (B/C) width. Latex free.

Black/Blue	Chocolate/Orange	Khaki/Red	Grey/Yellow	Women's Shoe Size
0815-89-464	0815-89-530	0815-89-605	0815-89-670	5
0815-89-472	0815-89-548	0815-89-613	0815-89-688	6
0815-89-480	0815-89-555	0815-89-621	0815-89-696	7
0815-89-498	0815-89-563	0815-89-639	0815-89-704	8
0815-89-506	0815-89-571	0815-89-647	0815-89-712	9
0815-89-514	0815-89-589	0815-89-654	0815-89-720	10
0815-89-522	0815-89-597	0815-89-662	0815-89-738	11



## Vionic Mojave. ⚠

Performance women's sandal with flexible leather uppers and two-tone webbing incorporates an elasticized toe ring plus hook-and-loop adjustable strap for a customized fit. Available in whole sizes only in medium/wide (B/C) width. Latex free.

Charcoal/Berry	Grey/Yellow	Women's Shoe Size
0815-96-162	0816-15-079	6
0815-96-170	0816-15-087	7
0815-96-188	0816-15-095	8
0815-96-204	0816-15-103	9
0815-96-212	0816-15-111	10



## Vionic Muir. ⚠

Performance women's back-strap sandal with flexible leather uppers and two-tone webbing incorporates hook-and-loop adjustable strap in both the forefoot and midfoot. The design offers a customized fit. Available in whole sizes only in medium/wide (B/C) width. Latex free.

Charcoal/Blue	Taupe	Women's Shoe Size
0815-96-220	0815-96-279	6
0815-96-238	0815-96-287	7
0815-96-246	0815-96-295	8
0815-96-253	0815-96-303	9
0815-96-261	0815-96-311	10



## Vionic Atlas. ⚠

Athletic style and support have evolved into the Atlas, the ideal toe post sandal for recovery after sports or a long day on foot. Features built-in orthotic footbed and dual density midsole for support and comfort. Available in whole sizes only in medium (B) width. Latex free.

Black/Silver	Grey/Coral	Women's Shoe Size
0815-95-966	0815-96-014	6
0815-95-974	0815-96-022	7
0815-95-982	0815-96-030	8
0815-95-990	0815-96-048	9
0815-96-006	0815-96-055	10



## Vionic Sprint. ⚠

Sprint to recovery with ease in this slip-on, after-sports sandal. With a pop of color for shelf appeal and femininity, plus medial side gore for individualized comfort and a triple-dense midsole, it rewards the foot with support and stability at every finish line. Available in whole sizes only in medium (B) width. Latex free.

Black/Silver	White/Citron	Women's Shoe Size
0815-96-063	0815-96-113	6
0815-96-071	0815-96-121	7
0815-96-089	0815-96-139	8
0815-96-097	0815-96-147	9
0815-96-105	0815-96-154	10

⚠ This product contains natural rubber latex or dry natural rubber.

Half sizes for men should go up a size (i.e. 10½ goes to 11)  
Half sizes for women should go down a size (i.e. 7½ goes to 7)



# Orthaheel Footwear

Signature Orthaheel Comfort.

**orthaheel**  
well-being with every step



## Orthaheel® Tide. ⚠

Fun and comfortable, the favorite Tide flip flop sandal makes it easy to support your feet and add color to your wardrobe as the weather warms up. Available in Women's whole sizes in three fashionable colors! Latex free.

Black	Bronze Metallic	Pewter Metallic	Women's Shoe Size
0815-96-659	0815-96-725	0815-96-790	5
0815-96-667	0815-96-733	0815-96-808	6
0815-96-675	0815-96-741	0815-96-816	7
0815-96-683	0815-96-758	0815-96-824	8
0815-96-691	0815-96-766	0815-96-832	9
0815-96-709	0815-96-774	0815-96-840	10
0815-96-717	0815-96-782	0815-96-857	11



## Orthaheel® Ryder. ⚠

A classic flip flop sandal with built-in orthotic support and added adjustability. Ryder also features a padded jersey-lined upper and a soft nylon webbing toe post for even greater comfort. Available in Men's whole sizes in two classic colors. Latex free.

Black/Grey	Chocolate/Tan	Men's Shoe Size
0815-96-865	0815-96-931	7
0815-96-873	0815-96-949	8
0815-96-881	0815-96-956	9
0815-96-899	0815-96-964	10
0815-96-907	0815-96-972	11
0815-96-915	0815-96-980	12
0815-96-923	0815-96-998	13



## Orthaheel® Unisex Wave Flip Flops. ⚠

The popular, iconic flip flop sandal from Vasyli features the built-in biomechanical Orthaheel footbed. The lightweight "heel" supports your patients' feet while helping to realign their lower legs and improve posture. Durable rubber outsole with wave-patterned tread helps improve traction. Latex free.

Unisex Wave Flip Flop				Men's Shoe Size	Women's Shoe Size
Black	Chocolate	Khaki	White/Blue		
5647-09	0815-53-668	5646-99	5647-19	5	6
5647-10	0815-53-676	5647-00	5647-20	6	7
5647-11	0815-53-684	5647-01	5647-21	7	8
5647-12	0815-53-692	5647-02	5647-22	8	9
5647-13	0815-53-700	5647-03	5647-23	9	10
5647-14	0815-53-718	5647-04	5647-24	10	11
5647-15	0815-53-726	5647-05	5647-25	11	12
5647-16	0815-53-734	5647-06	—	12	13
5647-17	0815-53-742	5647-07	—	13	14

Half sizes for men should go up a size (i.e. 10½ goes to 11)  
Half sizes for women should go down a size (i.e. 7½ goes to 7)

⚠ This product contains natural rubber latex or dry natural rubber.

Signature Orthaheel Comfort.

**orthaheel**  
well-being with every step



## Orthaheel® Relax Women's Slippers ⚠.

Soft and comfy women's slipper with terry uppers and hook-and-loop adjustable straps. This style can be worn year-round for much-needed support around the house. Fabric-covered TPR outsole with wave-pattern tread designed for indoor use. Built into the slipper is our award-winning biomechanical orthotic footbed which supports your feet and helps realign the lower legs—improving posture and balance. EcoFresh® antibacterial, odor-reducing footbed and linings. Latex free.

Black	Tan	Zebra	Women's Shoe Size
0815-80-588	0815-48-833	0816-17455	5.5-6
0815-80-596	0815-48-841	0816-17471	6.5-7
0815-80-604	0815-48-858	0816-17489	7.5-8
0815-80-612	0815-48-866	0816-17497	8.5-9
0815-80-620	0815-48-874	0816-17505	9.5-10

⚠ This product contains natural rubber latex or dry natural rubber.



## Orthaheel® Relief ¾ Length.

EVA orthotic with Tri-Planar Motion Control Technology. Cushioned shock dot in the heel for added shock absorption. Cut-away on the outside edge making it easy to fit into most mens and womens flat or low heeled shoes, or shoes without removable insole. EcoFresh® anti-bacterial odor absorbing cover. Contoured around the heel and arch areas to achieve 100% foot contact.



## Orthaheel® Relief Full Length.

EVA orthotic with Tri-Planar Motion Control Technology. Cushioned shock dot in the heel for added shock absorption. Can be trimmed in forefoot if necessary. 4-degree rear foot wedge to provide support and control which helps prevent over-pronation. EcoFresh® anti-bacterial odor absorbing cover. Contoured around the heel and arch areas to achieve 100% foot contact.



## Orthaheel® Active Full Length.

EVA orthotic with re-enforced hardened plastic. Cushioned shock dot in the heel for added shock absorption. Can be trimmed in forefoot if necessary. Designed to provide greater control in faster paced activities such as running and fast walking. 4-degree rear foot wedge to provide support and control which helps prevent over-pronation. EcoFresh® anti-bacterial odor absorbing cover. Contoured around the heel and arch areas to achieve 100% foot contact.

Relief 3/4 Length	Relief Full Length	Active Full Length	Size	Men's Shoe Size	Womens Shoe Size
0815-93-300	0815-93-250	0815-93-201	XS	3.5-5	4.5-6
0815-93-318	0815-93-268	0815-93-219	S	5.5-7	6.5-8
0815-93-326	0815-93-276	0815-93-227	M	7.5-9	8.5-10
0815-93-334	0815-93-284	0815-93-235	L	9.5-11	10.5-12
0815-93-342	0815-93-292	0815-93-243	XL	11.5-13	N/A



## Orthaheel® Gel Heel Device.

Made from soft PU Gel. Designed to fit almost all types of footwear. Great shock-absorbing qualities with EcoFresh® ant-bacterial deodorizer top cover. Softer rear-foot wedge provides low-level control of over-pronation

Gel Heel Device	Size	Men's Shoe Size	Women's Shoe Size
0815-93-359	Small	5-7.5	6-8.5
0815-93-367	Medium	8-10.5	9-11.5
0815-93-375	Large	11-13	12-14.5

Half sizes for men should go up a size (i.e. 10½ goes to 11)  
Half sizes for women should go down a size (i.e. 7½ goes to 7)



Coming in Fall

Look for our new line of footwear



**VASYLI®**

*Performance through innovation*

Contact your Patterson Medical Sales Representative for more information or visit [www.vasylimedical.com](http://www.vasylimedical.com) to view upcoming seminar dates and locations.






**Vasyli Medical Orthotic Patient Education Brochure**  
Package of 25  
0815-28-348



**Vasyli 4-Fold Family Product Flyer**  
Package of 25  
0816-15-145



**Vionic Sandal Flyer**  
Package of 25  
0816-14-148



**Orthaheel Sandal Flyer**  
Package of 25  
0816-13-967









**Orthaheel Insole Flyer**  
Package of 25  
0816-02-770

**Orthaheel Sandal Poster**  
(22" x 28")  
0816-14-122


**Grid Display**  
Versatile, can be used as a floor stand or a wall mount display.  
\*0816-12-563  
\* Item number reflects our 2013 Spring Kiosk. See your Patterson Medical sales representative for new footwear selections.



# Vasyli Medical Product Matrix

CATEGORIES	Armstrong Sensitive Feet	Dananberg 1st Ray Function	Hoke Supination	McConnell High Heels	McPoil Tissue Stress	Prior Sports
     						
Red Box indicates application						
CONDITION						
Back Pain						
Knee Pain						
Foot/Heel Pain				Arch/Forefoot		
Diabetic						
USAGE						
Running						
Walking						
Day to Day						
Lateral Activity (Basketball etc.)						
After Activity/Early Morning						
SHOE TYPE						
Athletic						Low Profile
Men's Casual	Depth Essential					
Women's Casual	Depth Essential					
Men/Ladies Casual (Narrow)						
Women's Dress (High Heel)						
Men's Dress						
Custom Fit Shoes						
Kids Shoes						
Football Boots		Deep Heel In Shoe	Deep Heel In Shoe			Low Profile
Work Boots						
Work Boots (High Heat Environment)						
SIZES AVAILABLE	XS, S, M, L, XL	XS, S, M, L, XL	XS,S,M,L,XL	XS,S,M,L,XL	XS,S,M,L,XL	XS, S, M, L, XL
PRODUCT AVAILABLE	NOW	NOW	September 2013	NOW	September 2013	NOW



CATEGORIES	Custom Range 3/4 & Full Length	Easy Fit	Extended Slimfit	Shock Absorber	Heel Pain Stabilizer	Orthaheel/Vionic Sandals
<div>Red Box indicates application</div> 						
CONDITION						
Back Pain						
Knee Pain						
Foot/Heel Pain						
Diabetic						
USAGE						
Running						
Walking						
Day to Day						
Lateral Activity (Basketball etc.)						
After Activity/Early Morning						
SHOE TYPE						
Athletic						
Men's Casual						
Women's Casual						
Men/Ladies Casual (Narrow)						
Women's Dress (High Heel)						
Men's Dress						
Custom Fit Shoes						
Kids Shoes						
Football Boots			Shallow Heel		Shallow Heel	
Work Boots	Full Length					
Work Boots (High Heat Environment)						
SIZES AVAILABLE	KS,KL,XS,S,M,L,XL,XXL	KL,XS,S,M,L,XL	XS,S,M,L	XS,S,M,L,XL	S,L	5-13 US
PRODUCT AVAILABLE	NOW	NOW	NOW	NOW	September 2013	NOW



# Partnership Performance Innovation

Our feet were designed to walk on natural, soft surfaces like soil and sand. When the feet strike soil or sand, these surfaces give way to accommodate the lower limb's natural 6 degree angle. Unfortunately, we give our feet unnatural hard, flat surfaces like pavements and floors. This forces our feet to roll over and our arches to flatten – a condition called excess pronation. Poor alignment of the feet and legs can cause wear and tear to other parts of the body, often disrupting normal knee function and hip alignment as well as increasing forces on muscles in the lower back.

**Vasyli products help restore natural foot function.**



Available to Rehab Clinicians exclusively through Patterson Medical

  
**PATTERSON**  
MEDICAL  
[www.pattersonmedical.com](http://www.pattersonmedical.com)

Call toll free:  
**1-800-323-5547**  
In Canada call 1-800-665-9200

Fax toll free:  
**1-800-547-4333**  
In Canada fax 1-877-566-8611

©2013 Patterson Medical Products, Inc.  
All rights reserved. Printed in USA.  
081587476LIT

*Performance through innovation*

**VASYLI®**